
多糖分析服务 / Glycan Analysis Services: 唾液酸定量分析 / Quantitative Sialic Acid Analysis

 Ludger

*Sialic acid
Reference
Panel*

Ludger
DMB (1,2-diamino
dihydroxy)
Cat # LT-DMB
Batch # B557

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Quantitative Sialic Acid Analysis and why it is important

Sialic acid analysis is a regulatory requirement laid out in the ICH Q6B guidelines for characterisation of biopharmaceuticals.

Sialic acids are negatively charged monosaccharides found on the non-reducing termini of glycans. They are important for the stability and 3D conformation of glycoproteins and are involved in many biological interactions. Sialic acids often have a pivotal functional impact: for example sialylation of the N-glycans on IgG increases anti-inflammatory activity; O-acetylated sialic acids can change ligand interactions and affect degradation (O-acetylated sialic acids including Neu5,9Ac₂ are found on EPO); and the presence of sialic acids also increases the serum half life of glycoproteins by preventing uptake by the asialoglycoprotein-receptor located on liver cells.

A diverse range of sialic acids are found in nature, but the two major sialic acids found on N-glycans and O-glycans in biopharmaceuticals are N-acetyl-neuraminic acid (Neu5Ac, or NANA) and N-glycolyl-neuraminic acid (Neu5Gc, or NAGA). Humans cannot synthesise Neu5Gc and its presence on a drug can lead to immune reactions such as chronic inflammation. Anti-Neu5Gc antibodies have been detected in normal human sera, and can neutralize any Neu5Gc containing biopharmaceutical, thus lowering the drug's efficacy. It is important to be aware that the choice of cell line can greatly influence the type of sialic acids present on a biopharmaceutical, for instance a large proportion of the sialic acids on mouse IgG are often Neu5Gc.

It is therefore imperative, for drug safety and efficacy, to monitor both the level and types of sialic acids during all stages of the product life cycle as well as QC batch to batch consistency.

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Ludger's Method – using our LT-KDMB-A1 kit

Sialic acids are released
from the glycoprotein by mild acid hydrolysis
(2 hours incubation at 80°C with 2M acetic acid)

Samples are labelled with DMB
(reacts with the α -keto carboxylic acid which
is not present on neutral monosaccharides)


Samples are run on RP-LC

Data analysis.

The absolute amounts of Neu5Ac & Neu5Gc are calculated by reference to standard curves from quantitative Neu5Ac & Neu5Gc standards. The presence of O-acetylated sialic acids is determined by comparison to a Neu5,9Ac₂ standard and to a reference panel containing Neu5Ac, Neu5Gc, Neu5,7Ac₂, Neu5,Gc9Ac, Neu5,9Ac₂ and Neu5,7,(8),9Ac₃Gc. The relative proportions of these different sialic acids are determined from the peak areas.

唾液酸定量分析: EPO / Quantitative Sialic Acid Analysis: EPO

- 作为nmol/mg蛋白定量的唾液酸(Neu5Ac和Neu5Gc) / Quantitation of sialic acids (Neu5Ac & Neu5Gc) as nmol/mg protein
- 相对定量的O-乙酰化的Neu5,9Ac₂ / Relative quantitation of the O-acetylated Neu5,9Ac₂

参照标准曲线测定nmol/mg蛋白中Neu5Gc、Neu5Ac的含量 / Neu5Gc & Neu5Ac quantified in nmol/mg protein by reference to standard curves


	nmol/mg protein	
Sialic Acid	Average	%CV
Neu5Gc	3.61	5.11
Neu5Ac	270.32	5.31

Neu5Gc、Neu5Ac、Neu5,9Ac₂的相对比例按%峰面积计算。 / Relative proportions of Neu5Gc, Neu5Ac & Neu5,9Ac₂ calculated from % peak areas.

	Relative % Area	
Sialic Acid	Average	%CV
Neu5Gc	0.99	0.24
Neu5,Ac	93.91	0.02
Neu5,9Ac ₂	5.09	0.39

唾液酸定量分析: IgG / Quantitative Sialic Acid Analysis: IgG

- 参照标准曲线测定唾液酸(Neu5Ac和Neu5Gc)作为nmol/mg蛋白的含量 / Quantitation of sialic acids (Neu5Ac & Neu5Gc) as nmol/mg protein by reference to standard curves

Mouse IgG	nmol/mg protein	
Sialic Acid	Average	%CV
Neu5Gc	1.65	6.06
Neu5Ac	0.53	5.41

Human IgG	nmol/mg protein	
Sialic Acid	Average	%CV
Neu5Gc	0.00	-
Neu5Ac	8.90	2.81